


A close-up photograph of a hand holding a ureterorenoscope. The device is a complex piece of medical equipment with a black handle, a clear tube, and a white sheath. A green tube is inserted into the device, and a blue tube is also visible. The background is a blue surgical drape.

# URETERORENOSCOPES

# URETERORENOSCOPE - DIRECT VIEW



PRODUCT CODE	DIAMETER	LENGTH	LENS ANGLE
DS-30-0645-00	6 mm / 7.5 fr	315 mm	5°
DS-30-0631-00			12°
DS-30-0646-00	6 mm / 7.5 fr	430 mm	5°
DS-30-0626-00			12°
DS-30-0633-00	8 mm / 9.8 fr	315 mm	12°
DS-30-0632-00		430 mm	



# URETERORENOSCOPE - OFFSET EYEVIEW





PRODUCT CODE	DIAMETER	LENGTH	LENS ANGLE
DS-30-0684-00	4.5 mm / 6.5 fr	315 mm	5°
DS-30-0639-00	6 mm / 7.5 fr	315 mm	5°
DS-30-0558-00			12°
DS-30-0685-00	4.5 mm / 6.5 fr	430 mm	5°
DS-30-0640-00	6 mm / 7.5 fr	430 mm	5°
DS-30-0514-00			12°
DS-30-0069-00	8 mm / 9.8 fr	315 mm	12°
DS-30-0070-00		430 mm	12°
DS-30-0404-00	8.5 mm / 11.5 fr	315 mm	12°
DS-30-0405-00			12°

# URETERORENOSCOPE - TILTED EYEPIECE



PRODUCT CODE	DIAMETER	LENGTH	LENS ANGLE
DS-30-0682-00	4.5 mm / 6.5 fr	315 mm	5°
DS-30-0641-00	6 mm / 7.5 fr	315 mm	5°
DS-30-0552-00			12°
DS-30-0683-00	4.5 mm / 6.5 fr	430 mm	5°
DS-30-0642-00	6 mm / 7.5 fr	430 mm	5°
DS-30-0528-00			12°
DS-30-0395-00	8 mm / 9.8 fr	315 mm	12°
DS-30-0403-00		430 mm	
DS-30-0406-00	8.5 mm / 11.5 fr	315 mm	12°
DS-30-0407-00			



**LENS ANGLE**