




**EAR, NOSE, THROAT
SURGICAL TELESCOPES**

OTOSCOPE – SCOPES RANGE



PRODUCT CODE	DIAMETER	LENGTH	LENS ANGLE
DS-30-0041-00	2.7 mm	50 mm	0°
DS-30-0042-00			30°
DS-30-0043-00	4 mm	50 mm	0°
DS-30-0044-00			30°

SINUSCOPE – SCOPES RANGE



PRODUCT CODE	DIAMETER	LENGTH	LENS ANGLE
DS-30-0029-00	2.7 mm	110 mm	0°
DS-30-0030-00			30°
DS-30-0031-00			70°

PRODUCT CODE	DIAMETER	LENGTH	LENS ANGLE
DS-30-0142-00	2.7 mm	140 mm	0°
DS-30-0143-00			30°
DS-30-0032-00	2.7 mm	175 mm	0°
DS-30-0033-00			30°
DS-30-0034-00			70°
DS-30-0035-00	2.7 mm	187.5 mm	0°
DS-30-0036-00			30°
DS-30-0037-00			70°
DS-30-0038-00	4 mm	175 mm	0°
DS-30-0039-00			30°
DS-30-0311-00			45°
DS-30-0040-00			70°

LARYNGOSCOPE - RANGE



PRODUCT CODE	DIAMETER	LENGTH	LENS ANGLE
DS-30-0048-00	10 mm	190 mm	70°
DS-30-0049-00			90°

BRONCHOSCOPE- RANGE

PRODUCT CODE	DIAMETER	LENGTH	LENS ANGLE
DS-30-0050-00	5.5 mm	489 mm	0°
DS-30-0051-00			30°

KEY BENEFITS

- 1- Improves visibility and allows more precise and safe surgery
- 2- Realizes rich color reproduction and provides suitable colors for each clinical discipline.
- 3- Enables easier determination of tissue boundaries
- 4- Better visualization of blood vessels and lesions due to high range of red color reproduction
- 5- Allows surgeons to observe fine patterns and structures of tissues in the body in extreme detail even when enlarged.